

Physik


Chemie · Biologie

Technik

LEYBOLD

LEYBOLD DIDACTIC GMBH

4/94

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Polarization filters are used to produce linearly polarized light and to qualitatively investigate polarization phenomena.

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1 Safety note

If the polarization filter is placed directly in front of the light source, use a heat protection filter to prevent damage to the plastic foil from overheating.

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2 Description, technical data

The filters consist of dichroic plastic foils (40 mm dia.) embedded between glass plates. The dichroism is realized by the aligned, dyed filamentary molecules of the foils. The light waves with an electric vector parallel to the filamentary molecules are almost completely absorbed by the dye. An electric vector perpendicular to the filaments is only slightly attenuated. The degree of polarization is more than 99%, and within the visible range it is virtually independent of the wavelength.

The polarization filters are fitted in holders (13 cm dia.) so that they can rotate. The rotating part of the holder is equipped with a scale ($0^\circ \dots \pm 90^\circ$ with 5° scale divisions). The filters are fitted in such a manner that the electric vector of the transmitted light oscillates in the direction of the indicator.

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3 Use

The holders of the filters are of the same size as the lens holders, and therefore the filters can easily be used on the small optical bench (460 42). Generally the first filter is used as the polarizer and the second filter as the analyzer. They can also be used in combination with analyzers (e.g. mirror analyzer 472 79).

In view of the fact that the polarizing effect on rays diverging considerably from the optical axis is by nature weaker, the filters should be used in a parallel light beam.

Furthermore, care should be taken that the filters are not damaged by overheating. For a filter mounted immediately in front of the lamp, the risk of damage is particularly great in the case of longer experiments.

Recommended heat protection filters:

Picture slider with built-in heat protection filter (450 66) for halogen lamp 12 V/50 W/100 W (450 64) or plate glass cell (e.g. 477 20) filled with a weak solution of copper sulfate.